


**Shoulder Subluxation**


Shoulder Subluxation often results from the weight of an injured or flaccid arm applying direct mechanical stretch to the joint capsule of the shoulder. The weight and resultant stretch partial dislocate the joint and create decreased strength and impaired upper extremity function and coordination. Shoulder subluxation has been implicated as a causative factor for RSD. NMES has provided some moderate success in the prevention and treatment of subluxation. Research shows that, not only can NMES help reduce the subluxation of the shoulder joint, but it enhances motor recovery of the shoulder and also helps reduce shoulder pain.

Suggested Electrode Placements**2-pad placement**

Place electrode pads on the infraspinatus/teres minor and supraspinatus muscles of the shoulder

**4-pad placement**

Place electrode pads on the infraspinatus/teres minor, supraspinatus posterior deltoid, and anterior deltoid of the shoulder

**Suggested Settings**

Pulse rate: 50-65 Hz

Pulse width/duration: 250-300 μ s

On/Off times: 1:3 ratio (start with 5 sec ON, 15 sec OFF)

Ramp time: 1 sec

Amplitude: Maximally tolerated contraction

Treatment Time: 15-20 minutes

Sessions per day: 2-3 sessions per day (or as needed)

